

Introduction

Across all healthcare settings, preventing pressure injuries and being proactive by implementing early interventions has always been the gold standard of care. With recent Covid-19 isolation precautions and limited staff interaction, the rapid development of pressure injuries, along with falls and resident comfort is an ongoing challenge for many high-risk residents in long-term care facilities.

Five long-term care facilities participated in a 4-month product trial with a reactive air overlay* and documented the outcomes.

Objective

Prevent pressure injury development, avoid falls, promote comfort and provide pressure redistribution to high-risk residents in long-term care by implementing a reactive air overlay.

Method

Data was reviewed for high-risk residents according to their risk assessment. Skin assessments were performed by the wound team. Residents with current pressure injuries were excluded. We selected residents that had at least 1 of the following:

- High-risk according to their risk assessment
- Decreased oral intake
- Weight loss
- Isolation rooms
- Recent decline in health
- Any Covid-19 residents

An easy visual education tool was developed to determine the appropriate residents. 42 reactive air overlays were initially delivered to the 5 facilities. Training and in-service were provided. Residents were monitored weekly for any skin impairment and/or pressure injuries. Performed hand checks under the buttocks area to ensure proper air inflation.

* WAFFLE® Overlay is a registered trademark of EHOB, Inc. Indianapolis, IN

Challenges

Covid-19

- Resident’s health and mental well-being were declining quickly
- Staff were hesitant and reluctant to enter resident’s rooms
- Overall, resident/staff contact decreased

Staff Challenges

- Staffing greatly decreased for various reasons
- Shortage of Wound Coordinators and Wound Nurses

Remote Staff Training

- CDC guidelines mandated virtual training
- More challenges with virtual training vs. hands on training
- Staff was provided with 1-page educational sheet on how to use and apply the reactive air overlay

Study Obstacles

- Staffing challenges in the wound care department
- Due to turnover & limited staff, proper hand checks & air inflation was not performed routinely
- Overlay was not transferred consistently with residents as they moved from one room to another
- Product training could not keep up with staff turnover

Results – Prevented Pressure Injuries • No Falls • Enhanced Patient Comfort

Joliet (16 on Overlay)

- Zero skin breakdown in 15 residents.
- 1 resident on hospice developed a Stage 2 heel ulcer History included Adult Failure to Thrive, Coronary Artery Disease, Anemia and Dementia Resident came into the facility and was placed on a reactive air overlay and heel protectors. Despite all the proactive interventions, along with frequent repositioning, the resident developed a Stage 2
- Zero falls

South Shore (35 on Overlay)

- Zero skin breakdown
- Zero falls

Morgan Park (35 on Overlay)

- Zero skin breakdown
- Zero falls
- Residents reported liking the overlay

Bronzeville (36 on Overlay)

- Zero skin breakdown
- Zero falls

Aria (36 on Overlay)

- Zero skin breakdown in 35 residents
- 1 resident develop a Stage 2 to the right buttock area. He was admitted with MASD to the bilateral buttock/sacral area. Incontinent of bowel and bladder. Multiple co-morbidities which included history of Covid-19, CVA, protein-calorie malnutrition, anemia and developed a stage 2 just 1 week after being admitted to Aria.
- Zero falls

Additional Study Feedback

- Residents appreciated the comfort of the reactive air overlay
- Non-powered, adding an advantage of no tubing on the floor
- Decreased the usage/rentals of low air Loss mattress – financial benefit
- Prevented shear and friction with transfers and transportation (i.e. dialysis, outpatient appointments, etc.)

Conclusions

- **Out of 158 residents, only 2 pressure injuries developed**
- **Zero falls were reported**
- **Patient comfort increased**
- **Excellent staff & resident satisfaction**
- **Cost savings achieved using this overlay vs. the rental of a low air loss mattress**

As a result of this trial, we have incorporated the use of this overlay across all our facilities. Additional follow-up is recommended to determine long term benefits.

Discussion

Recent Covid-19 challenges were detrimental to resident care. Implementing early interventions with the reactive air overlays resulted in only two Stage 2 pressure injuries, suggesting that early interventions for high-risk residents can decrease pressure injuries and improve quality of life.

