



Indications for Use:

- For lifting the heel off of the bed for patients who are at risk of heel pressure injury.
- For assisting in preventing and managing heel pressure injuries in patients who have existing or history of heel pressure injuries.
- For assisting in preventing skin injury to the foot for patients at risk.
- For assisting in preventing hip rotation with use of optional wedge.

Contraindications

- Do not use for ambulation or chair transfers.

Interventions

- Remove each shift and, as needed, for skin assessments.
- Readjust to ensure proper placement after repositioning the patient.


Instructions for Use




Lower the foot into the product, ensuring the heel is centered over the heel opening.




Close the foot strap, then the leg straps.



Confirm proper fit by inserting your fingertips into the product along the side of the foot. Your fingers should slide into the opening without difficulty. If they do not slide easily, loosen the foot strap and recheck.



Perform a hand check by sliding hand, palm up, under the heel to verify that the heel is offloaded. A handcheck should be performed every two hours and with patient repositioning.



The Open Foot Gate can be used for skin assessments, for verifying heel offloading, or for accommodating a wider foot.



Optional Anti-Rotation Wedge may be placed on the side of the product to help prevent hip and leg rotation.

*Tubing ports may be used to route tubing away from the skin.

Cleaning Instructions

- Cleanable with common hospital disinfectant wipes that do not contain bleach.

Instructional Video* ►



For complete instructions, refer to the package insert.

**For information on sizing and tubing port usage, please watch video.*