



Indications for Use

The TruVue® Heel Protector is a single patient use heel offloading device. It helps to prevent pressure injuries, while reducing the risk for plantarflexion contractures and maintaining the ankle in a neutral position (optional wedge recommended).


Contraindications

Do not use for ambulation or chair transfers.

Interventions

Remove each shift and, as needed, for skin assessments. Readjust to ensure proper placement after repositioning the patient.

Instructions for Use




1. Align Foot

Align the patient's ankle bone with the top of the heel opening.



2. Assemble

Fold the sides of the heel protector up to assemble the Open Foot Gate.



3. Secure Straps

Secure the Flexible Straps across the leg and foot. Adjust the Anti-Foot Drop Straps as needed to help keep the foot at a 90 degree angle.




4. Attach Optional Wedge

Attach the Anti-Rotation Wedge to the side of the heel protector to help prevent hip and leg rotation.



5. Verify Offloading

Place your hand under the heel to ensure that the heel is offloaded.



6. Skin Assessments

The Open Foot Gate can be used for skin assessments, to verify the heel is offloaded, or for patient comfort.

Cleaning Instructions

The TruVue Heel Protector can be cleaned with the standard hospital cleaning wipes if soiled. A complete list of approved cleaning products is available upon request.

Instructional Video*



For complete instructions, refer to the package insert.

*For information on sizing and tubing port usage, please watch video.