

# WAFFLE®

## Cushions



### Pre-Inflated Cushions

WAFFLE® Cushions are shipped pre-inflated to their ideal volume and should only appear to be 60% full. This allows for the patient to immerse into the product, while adequately redistributing pressure over bony prominences.



Place the pre-inflated cushion with the label down and to the back of the seating surface. If the cushion has a valve, it should also be facing down.

### Indications for Use

- Prevention and management of pressure injuries
- Patient comfort
- Use across the continuum of care

### Precautions

- If used in a chair, position shifts should be made every 15–30 minutes or per individualized schedule
- Avoid direct contact between skin and product by using a pillowcase, sheet, or breathable under-pad
- Single patient use
- Weight Limit:
  - Small: 100 lbs.
  - Medium: 300 lbs.
  - Large: 400 lbs.
  - X-Large: 700 lbs.

### Additional Uses



Behind Head



Under Elbow



Between Knees



Additional Support

### Cleaning Instructions

The WAFFLE® Cushion is single patient use but can be cleaned with standard hospital cleaning wipes if soiled. A complete list of approved cleaning products is available upon request.

**For complete instructions, refer to product package insert.**



Scan the QR code to view an instructional video.