

WAFFLE™

Cushion



Indications for Use

- Prevention and management of pressure injuries
- Patient comfort
- Use across the continuum of care

Precautions

- If used in a chair, position shifts should be made every 15–30 minutes or per individualized schedule
- Avoid direct contact between skin and product by using a pillowcase, sheet, or breathable under-pad
- Single patient use
- Weight Limit: 300 lbs.

Set-Up

Step 1: Inflate



Using the WAFFLE® L.A.D. Pump, insert the tip of the pump into the valve. Refer to the inflation chart on the side of the pump, and inflate your cushion using the recommended level of strokes.

Step 2: Check



Once the cushion appears about 60% full, check for proper inflation. If you can easily roll one side of the product past the first set of holes, but not to the second, your cushion is ready to use.

Step 3: Utilize



Place product on seating surface, valve down and towards the back of the surface.

WAFFLE® Cushion Additional Uses for Protection

Behind Head



Under Elbow



Between Knees



Additional Support



Cleaning Instructions

The WAFFLE Cushion is single patient use but can be cleaned with standard hospital cleaning wipes if soiled. A complete list of approved cleaning products is available upon request.

For complete instructions, refer to product package insert.



Scan the QR code to view an instructional video.



250 N Belmont Avenue • Indianapolis, IN 46222 • 800.899.5553 • EHO.com

WAFFLE is a trademark of EHO, Inc. © 2020 EHO, Inc. • MKG002427_REVo1_WC_Valve_Product_Education_Flyer_06-20_V2